

ATTACHMENT

Animal Test Kit PetAge360

A personalized comprehensive Health Check for companion animals, collected with an oral swab without blood draw or anesthesia

3 tests

Test items

No blood · No anesthesia

Noninvasive oral swab collection

6 steps

Request form → courier shipment

Room temperature

Storage & shipping

Product	PetAge360® Personalized Comprehensive Health Check for Companion Animals
Made & sold by	DNAGENICS
Test items	Telomere (aging rate) · Oral microbiome · Tumor genes
Sample	Collected by rubbing the inner gums with an oral swab (no blood draw or anesthesia)

01 Kit contents



Kit contents

Instruction leaflet · DNA EVER MICROBI box · Swab/preservative tube



Sample collection tube

Oral swab and preservative tube

02 Testing procedure (per the instruction leaflet)

STEP 1

Fill out the test request form.

STEP 4

Attach barcode stickers to the preservative tube and the request form.

STEP 2

Rub the swab along the gums as far as the inside of the canine teeth to collect the sample (saliva).

STEP 5

Place the request form and the preservative tube in the box provided.

STEP 3

Put the swab into the preservative tube and snap it off at the middle joint.

STEP 6

Ship by courier to the address on the side of the kit.

03 Precautions for sample collection

Checks before and after collection

- Do not let the animal eat for 30 minutes before collection.
- Do not let the animal drink water for 10 minutes before collection.
- Ship immediately after collection and store at room temperature.
- Make sure the preservative tube is tightly closed.
- Take care that the preservative solution is not swallowed.